

PTM/HT cartridge filters are made with Polytetrafluoroethylene (PTFE) membrane and high temperature polypropylene for cartridge hardware and membrane support. The filters are validated for sterilizing in elevated temperature gas and non-aqueous liquid applications (up to 221°F). Pore sizes range from 0.10 to 1.0 μm .


The hydrophobic PTM/HT filters resist wetting by airborne water droplets, making them ideal for air and gas applications. The broad chemical compatibility of the PTM filters makes them well suited for aggressive solvents and other non-aqueous liquids. Each cartridge module is individually tested using the water intrusion method before it is released from manufacture.

Critical Process provides unrivaled delivery times, technical consulting before purchasing, and very competitively priced high-performance products. Our comprehensive testing & analysis and validation services support your team whenever they need it. Your process experts partnering with our filtration experts is how we deliver your company's solution right the first time.

Sterilizing Filters Tank Vent & Process Gas



CARTRIDGES – Nominal Dimensions Length: 5 to 40 in. (12.7 to 101.6 cm) Outside Diameter: 2.75 in. (7.0 cm)



PTM/HT sterilizing filters are recommended for:

- Compressed Air
- Pressurized Gases
- Fermentation Air
- Solvents

Maximum Operating Parameters

Operating Temperature (air/gas)221 °F (105 °C)Forward Differential Pressure80 psid at 68 °F (5.52 bard at 20 °C) (Liquid and Gas)		CARTRIDGES
Forward Differential Pressure 80 psid at 68 °F (5.52 bard at 20 °C) (Liquid and Gas)	Operating Temperature (air/gas)	221 °F (105 °C)
	Forward Differential Pressure	80 psid at 68 °F (5.52 bard at 20 °C) (Liquid and Gas)
Reverse Differential Pressure 50 psid at 68 °F (3.45 bard at 20 °C)	Reverse Differential Pressure	50 psid at 68 °F (3.45 bard at 20 °C)

Sanitization & Sterilization

Filtered Hot Water*	90 °C (194 °F), 30 minutes, multiple cycles, max 3 psid forward flow
Inline Steam*	275 °F (135 °C), 30 min, 25+ cycles
Autoclave*	250 °F (121 °C), 30 min, 25+ cycles
Chemical Sanitization	Performed using industry standard concentrations of hydrogen peroxide, peracetic acid, sodium hypochlorite and other selected chemicals.

^{*}For all elevated temperature procedures above, a stainless-steel support ring is required.

Filtration Area (Nominal)

	CARTRIDGES				
Length	5"	10"	20"	30"	40"
	12.7cm	25.4cm	50.8cm	76.2cm	101.6cm
Area	3.8 ft ²	8.2 ft ²	16.4 ft ²	24.6 ft ²	32.8 ft ²
	0.36m ²	0.76m ²	1.52m ²	2.29m ²	3.05m ²

Integrity Testing

PORE SIZE	WATER INTRUSION TEST PRESSURE			E POINT MUM*
μm	PSIG	BARG	PSIG	BARG
0.10	35	2.41	22	1.52
0.22	35	2.41	18	1.24
0.45	N/A	N/A	9	0.62
1.0	N/A	N/A	6	0.41

WATER INTRUSION SPECIFICATIONS (mL/10 min)					
Length	5"	10"	20"	30"	40"
0.10μm	≤ 4.3	≤ 10	≤ 20	≤ 30	≤ 40
0.22μm	≤ 5.6	≤ 13	≤ 26	≤ 39	≤ 52

^{*} Bubble Point for membrane wetted with 60% IPA / 40% water solution.

Construction Materials

Filtration Media Polytetrafluoroethylene (PTFE) Membrane	
Media Support High Temperature Polypropylene	
End Caps, Center Core, Outer Support Cage	High Temperature Polypropylene
Sealing Method	Thermal Bonding
O-Rings/Gaskets Cartridges only Buna, Viton® (or FKM), EPDM, Sil FEP Encapsulated Silicone, FEP Encapsulated Viton (or FKM)	

Validation

PTM/HT filters are validated using test procedures that comply with ASTM F 838-15(ae1) protocols for the determination of bacterial retention in filters used for liquid filtration. The challenge level is a minimum of 10^7 organisms per \mbox{cm}^2 of filter media. CPF filters have > 7-log removal when challenged with the organisms listed below (0.10 μ m and 0.22 μ m meet the FDA definition of sterilizing grade filters).

0.10μm: Brevundimonas diminuta 0.22μm: Brevundimonas diminuta 0.45μm: Serratia marcescens

Validation Guides available upon request.

Endotoxins

The levels of bacterial endotoxins in aqueous extracts from PTM/HT filters are below current USP limits as specified for water for injection.

Extractables

PTM/HT filters typically exhibit low levels of non-volatile residues.

Toxicity Compliance

Materials used to construct the PTM/HT filters are nontoxic and meet the requirements for the MEM Elution Cytotoxicity Test and the requirements for Biological Reactivity Tests in the current version of the United States Pharmacopeia (USP) for Class VI - 121 °C Plastics.


Non-Fiber Releasing

PTM/HT filters comply with Title 21 CFR sections 210.3 (b)(6) and 211.72, for non-fiber releasing filters.


FDA Compliance

Materials meet the requirements listed by the FDA as appropriate for use in articles intended for repeated food contact as specified in Title 21 CFR sections 174.5, 177.1500, 177.1520, 177.1630, 177.2440, and 177.2600 as applicable.

Water Flow Rates for PTM/HT Cartridges by Pore Size



Air/Gas Flow Rates for PTM/HT Cartridges by Pore Size




Flow rates for Cartridge filters are per 10-inch length. The test fluid is water or air at ambient temperature.

PTM/HT Filters Ordering Information

Fill in the corresponding codes in the boxes below to build your Part Number.

To consult with one of our technical team members, request a quote or place an order: call (603) 880-4420 or contact us here.





One Chestnut Street Nashua, NH 03060 603.880.4420 FAX: 603.880.4536

CriticalProcess.com

The information contained herein is subject to change without notice. The Critical Process Filtration logo is a trademark of Critical Process Filtration, Inc. Viton is a trademark of DuPont Performance Elastomers L.L.C.

© 2020 Critical Process Filtration, Inc. • All Rights Reserved

Data Sheet PTMHTDS Rev B