


BCWPS cartridge and capsule filters are CPF's proprietary, high-capacity clarification and prefiltration solutions for particle retention in aqueous fluids. Their retention capacity exceeds that of traditional depth filters and is the better choice for fluids with heavy contamination. BCWPS filters deliver high flow and throughput across a wide pH range. The filters' low binding characteristics are an important feature when filtering fluids with high proteins and preservatives.

The BCWPS filters consist of a single layer of Polyethersulfone (PES) membrane. For very high particle loads, an optional prefilter made of the same membrane can be integrated. To achieve targeted results, you choose the pore sizes of each layer. Final layer pore size options range from 0.1 to 1.0 μ m, and prefilter pore size options range from 0.2 to 1.0 μ m

When multiple prefiltration steps are required in your process, the BCWPS dual layer option acts as two prefilters in one device, which reduces the amount of housings needed in your filtration train. CPF filter devices scale from laboratory to full production using identical materials to ensure consistent results. BCWPS capsule filters are available pre-sterilized.

Critical Process provides unrivaled delivery times, technical consulting before purchasing, and very competitively priced high-performance products. Our comprehensive testing & analysis and validation services support your team whenever they need it. Your process experts partnering with our filtration experts is how we deliver your company's solution right the first time.



BCWPS is the recommended media for clarification and prefiltration in:

- SVPs & LVPs
- Diagnostics
- Buffers
- WFI, Water Purification
- Vaccines
- Ophthalmics

Clarification & Prefiltration



CARTRIDGES – Nominal Dimensions Length: 5 to 40 in. (12.7 to 101.6 cm) Outside Diameter: 2.75 in. (7.0 cm)



CAPSULES – Nominal Dimensions Length: 2 to 30 in. (5.1 to 76.2 cm) Outside Diameter: 3.50 in. (8.9 cm)

Maximum Operating Parameters

| | CARTRIDGES | CAPSULES | |
|--------------------------------|---|--|--|
| iquid Operational Pressure | N/A | 80 psi at 68 °F (5.52 bard at 20 °C) | |
| Gases Operational Pressure | N/A 60 psi at 68 °F (4.14 bar at 20 °C) | | |
| Operating Temperature (water) | 180 °F at 30 psid (82 °C at 2.07 bard) 110 °F at 30 psid (43 °C at 2.07 bard) | | |
| Forward Differential Pressure | 80 psid at 68 °F (5.52 bard at 20 °C) | 8°F (5.52 bard at 20 °C) 80 psid at 68 °F (5.52 bard at 20 °C) | |
| Reverse Differential Pressure | 50 psid at 68 °F (3.45 bard at 20 °C) | 50 psid at 68 °F (3.45 bard at 20 °C) | |
| Recommended Changeout Pressure | 35 psid (2.41 bard) | 35 psid (2.41 bard) | |

Sanitization & Sterilization

| Filtered Hot Water* | 90 °C (194 °F), 30 minutes, multiple cycle max 3 psid forward flow | rs, N/A |
|-----------------------|--|---|
| Inline Steam* | 275 °F (135 °C), 30 min, 25+ cycles | N/A |
| Autoclave* | 250 °F (121 °C), 30 min, 25+ cycles | 250 °F (121 °C), 30 min, 25+ cycles |
| Chemical Sanitization | Performed using industry standard conce hypochlorite and other selected chemica | entrations of hydrogen peroxide, peracetic acid, sodium ls. |

*Cartridge Filters – For all elevated temperature procedures above, a stainless-steel support ring is required.

| | CAPSULES | CARTRIDGES AND CAPSULES | | | CARTRIDGES | |
|--------------|---------------------|-------------------------|---------------------|----------------------|----------------------|----------------------|
| Longth | 2″ | 5″ | 10" | 20" | 30" | 40" |
| Length | 5.08cm | 12.7cm | 25.4cm | 50.8cm | 76.2cm | 101.6cm |
| Area – | 1.1 ft ² | 3.1 ft ² | 6.7 ft ² | 13.4 ft ² | 20.1 ft ² | 26.8 ft ² |
| Single Layer | 0.11m ² | 0.29m ² | 0.62m ² | 1.24m ² | 1.86m ² | 2.48m ² |
| Area – | 0.9 ft ² | 2.5 ft ² | 5.3 ft ² | 10.6 ft ² | 15.9 ft ² | 21.2 ft ² |
| Dual Layer | 0.08m ² | 0.23m ² | 0.49m ² | 0.98m ² | 1.47m ² | 1.96m ² |

Filtration Area (Nominal)

Construction Materials

| Filtration Media | Single or Dual Layer High Capacity PES Membrane with Polyester Support | | |
|--|--|--|--|
| Media Support | Polypropylene | | |
| End Caps, Center Core, Outer Support Cage, Capsule Housing | Polypropylene | | |
| Sealing Method | Thermal Bonding | | |
| O-Rings/Gaskets Cartridges only | Buna, Viton [®] (or FKM), EPDM, Silicone, FEP Encapsulated Silicone, FEP Encapsulated Viton (or FKM) | | |

Endotoxins

The levels of bacterial endotoxins in aqueous extracts from BCWPS filters are below current USP limits as specified for water for injection.

Extractables

BCWPS filters typically exhibit low levels of non-volatile residues.

TOC and Conductivity

BCWPS filters conform with TOC standards of USP <643> and the water conductivity standards of USP <645> after an appropriate flush with purified water.

Toxicity Compliance


Materials used to construct BCWPS filters are non-toxic and meet the requirements for the MEM Elution Cytotoxicity Test and the requirements for Biological Reactivity Tests in the current version of the United States Pharmacopeia (USP) for Class VI - 121 °C Plastics.

Non-Fiber Releasing

BCWPS filters comply with Title 21 CFR sections 210.3 (b)(6) and 211.72, for non-fiber releasing filters.


FDA Compliance

Materials meet the requirements listed by the FDA as appropriate for use in articles intended for repeated food contact as specified in Title 21 CFR sections 174.5, 177.1500, 177.1520, 177.1630, 177.2440, and 177.2600 as applicable.



Flow Rates for BCWPS Cartridges by Pore Size

Flow rates for Cartridge filters are per 10-inch length. The test fluid is water at ambient temperature.




Flow Rates for BCWPS Capsules by Pore Size

Flow rates for Capsule filters are tested using a 2" capsule filter with 1" sanitary inlet and outlet ports. The test fluid is water at ambient temperature. Flow rates for larger capsules will scale with filtration area. Rates will vary based on end configuration of the capsule.

BCWPS Filters Ordering Information

Fill in the corresponding codes in the boxes below to build your Part Number.

To consult with one of our technical team members, request a quote or place an order: call (603) 880-4420 or <u>contact us here.</u>





One Chestnut Street Nashua, NH 03060 603.880.4420 FAX: 603.880.4536

The information contained herein is subject to change without notice. The Critical Process Filtration logo is a trademark of Critical Process Filtration, Inc. Viton is a trademark of DuPont Performance Elastomers L.L.C. © 2025 Critical Process Filtration, Inc. • All Rights Reserved

CriticalProcess.com