

FPS filters provide a wide selection of validated, single layer Polyethersulfone (PES) cartridge and capsule filters used for removing bacteria, mold and yeast in aqueous liquids. FPS filters have been designed to comply with all FDA requirements for the food industry. Pore sizes range from 0.03 to 1.2 μ m and the filter sizes scale from laboratory to full production using identical materials to ensure consistent results.

FPS filters have low protein binding characteristics making them a good choice for fermented beverage filtration. These filters deliver high throughput and handle cleaning and sanitization protocols well.


FPS filters are flushed to remove manufacturing debris and reduce extractables. Products are 100% integrity tested.

Critical Process provides unrivaled delivery times, technical consulting before purchasing, and very competitively priced highperformance products. Our comprehensive testing & analysis and validation services support your team whenever they need it. Your process experts partnering with our filtration experts is how we deliver your company's solution right the first time.

Bacteria/Yeast/Mold Removal



CARTRIDGES – Nominal Dimensions Length: 5 to 40 in. (12.7 to 101.6 cm) Outside Diameter: 2.75 in. (7.0 cm)



FPS filters are recommended for the filtration of:

- Wine
- Beer
- Clear Juices
- Aseptically Packaged Liquids
- Bottled Water
- Process Water
- Clean-in-Place Solutions



CAPSULES – Nominal Dimensions Length: 2 to 30 in. (5.1 to 76.2 cm) Outside Diameter: 3.50 in. (8.9 cm)

Maximum Operating Parameters

	CARTRIDGES	CAPSULES	
Liquid Operational Pressure	N/A	80 psi at 68 °F (5.52 bard at 20 °C)	
Gases Operational Pressure	N/A	60 psi at 68 °F (4.14 bar at 20 °C)	
Operating Temperature (water)	180 °F at 30 psid (82 °C at 2.07 bard)	110 °F at 30 psid (43 °C at 2.07 bard)	
Forward Differential Pressure	80 psid at 68 °F (5.52 bard at 20 °C) (Liquid and Gas)	Liquid - 80 psid at 68 °F (5.52 bard at 20 °C) Gas - 60 psi at 68 °F (4.14 bar at 20 °C)	
Reverse Differential Pressure	50 psid at 68 °F (3.45 bard at 20 °C)	50 psid at 68 °F (3.45 bard at 20 °C)	
Recommended Changeout Pressure	35 psid (2.41 bard)	35 psid (2.41 bard)	

Sanitization & Sterilization

Filtered Hot Water*	90 °C (194 °F), 30 minutes, multiple cycles, max 3 psid forward flow	N/A
Inline Steam*	275 °F (135 °C), 30 min, 25+ cycles	N/A
Autoclave*	250 °F (121 °C), 30 min, 25+ cycles	250 °F (121 °C), 30 min, 25+ cycles
Chemical Sanitization	Performed using industry standard co hypochlorite and other selected chem	ncentrations of hydrogen peroxide, peracetic acid, sodium nicals.

*Cartridge Filters – For all elevated temperature procedures above, a stainless-steel support ring is required.

Filtration Area (Nominal)

	CAPSULES	CARTRIDGES AND CAPSULES			CARTRIDGES	
Length	2″	5″	10"	20"	30"	40"
	5.08cm	12.7cm	25.4cm	50.8cm	76.2cm	101.6cm
Area	1.2 ft ²	3.4 ft ²	7.3 ft ²	14.6 ft ²	21.9 ft ²	29.2 ft ²
	0.11m ²	0.32m ²	0.68m ²	1.36m ²	2.04m ²	2.72m ²

Integrity Testing

PORE SIZE	DIFFUSION TEST PRESSURE*		BUBBLE POINT MINIMUM*	
μm	PSIG	BARG	PSIG	BARG
0.03	60	4.14	**	**
0.10	48	3.31	**	**
0.22	35	2.41	50	3.4
0.45	20	1.38	25	1.7
0.65	15	1.03	19	1.3
0.80	12	0.83	15	1.0
1.0	8	0.55	10	0.7
1.2	7	0.48	9	0.6

DIFFUSIC	ON SPECI	FICATION	S			
Length	2″	5″	10"	20"	30"	40"
mL/min	≤4.3	≤ 12.9	≤ 30	≤ 60	≤ 90	≤ 120

* For water wetted membrane

** Test pressure exceeds operational limits of capsule filters. Use the Diffusion Test method.

Construction Materials

Filtration Media	Single Layered Polyethersulfone (PES) Membrane
Media Support	Polypropylene
End Caps, Center Core, Outer Support Cage, Capsule Housing	Polypropylene
Sealing Method	Thermal Bonding
O-Rings/Gaskets Cartridges only	Buna, Viton® (or FKM), EPDM, Silicone, FEP Encapsulated Silicone, FEP Encapsulated Viton (or FKM)

Validation

FPS filters are validated using test procedures that comply with ASTM F 838-15(ae1) protocols for the determination of bacterial retention in filters used for liquid filtration. The filters are challenged with the organisms listed below.

0.03μm: Acholeplasma laidlawii 0.10μm: Brevundimonas diminuta 0.22μm: Brevundimonas diminuta 0.45μm: Serratia marcescens 0.65μm: Saccharomyces cerevisiae

Extractables


FPS filters typically exhibit low levels of non-volatile residues.

Non-Fiber Releasing

The FPS filters comply with Title 21 CFR sections 210.3 (b)(6) and 211.72, for non-fiber releasing filters.


FDA and EC Compliance

Materials meet the requirements listed by the FDA as appropriate for use in articles intended for repeated food contact as specified in Title 21 CFR sections 174.5, 177.1500, 177.1520, 177.1630, 177.2440, and 177.2600 as applicable. All materials used to make the filters are listed in European Commission Regulation EU/10/2011, Annex 1.



Flow Rates for FPS Cartridges by Pore Size

Flow rates for Cartridge filters are per 10-inch length. The test fluid is water at ambient temperature.




Flow Rates for FPS Capsules by Pore Size

Flow rates for Capsule filters are tested using a 2" capsule filter with 1" sanitary inlet and outlet ports. The test fluid is water at ambient temperature. Flow rates for larger capsules will scale with filtration area. Rates will vary based on end configuration of the capsule.

FPS Filters Ordering Information

Fill in the corresponding codes in the boxes below to build your Part Number.

To consult with one of our technical team members, request a quote or place an order: call (603) 880-4420 or <u>contact us here</u>





One Chestnut Street Nashua, NH 03060 603.880.4420 FAX: 603.880.4536

CriticalProcess.com

The information contained herein is subject to change without notice. The Critical Process Filtration logo is a trademark of Critical Process Filtration, Inc. Viton is a trademark of DuPont Performance Elastomers L.L.C. © 2025 Critical Process Filtration, Inc. • All Rights Reserved

Data Sheet FPSDS Rev D