



BCWPS Micro Capsule Filters


High Capacity PES Membrane

BCWPS Micro Capsule filters are constructed with a proprietary high capacity Polyethersulfone (PES) membrane and are available in single or dual layers. They are used for reducing particles in critical aqueous liquids and prefiltration in multiple applications. Pore sizes range from 0.1 to 1.0 μm and filter sizes scale from laboratory to full production using identical materials to ensure consistent results.

BCWPS Micro Capsule hydrophilic filters have a high contaminant holding capacity. They are utilized in clarifying and prefiltering products with high proteins and preservatives. The filter's low binding characteristics make them highly efficient, which is critical to protecting downstream filters from premature fouling and extending the life of final filters.

The BCWPS Micro Capsule filter delivers high flow and throughput across a wide pH range. BCWPS Micro Capsules are available pre-sterilized.

Critical Process provides unrivaled delivery times, technical consulting before purchasing, and very competitively priced high-performance products. Our comprehensive testing & analysis and validation services support your team whenever they need it. Your process experts partnering with our filtration experts is how we deliver your company's solution right the first time.



BCWPS Micro Capsule filters are the recommended media for clarification and prefiltration in:

- SVPs & LVPs
- Diagnostics
- Buffers
- WFI, Water Purification
- Vaccines
- Ophthalmics

Clarification & Prefiltration



MICRO CAPSULES – Nominal Dimensions

Body Length: 1.9 in. (4.8 cm)

Overall Length: 2.8 to 3.8 in. (7.1 to 9.7 cm)

Outside Diameter: 2.6 in. (6.6 cm)

Maximum Operating Parameters

MICRO CAPSULES	
Liquid Operational Pressure	80 psi at 68 °F (5.52 bard at 20 °C)
Gases Operational Pressure	60 psi at 68 °F (4.14 bar at 20 °C)
Operating Temperature (water)	110 °F at 30 psid (43 °C at 2.07 bard)
Forward Differential Pressure	80 psid at 68 °F (5.52 bard at 20 °C)
Reverse Differential Pressure	50 psid at 68 °F (3.45 bard at 20 °C)
Recommended Changeout Pressure	35 psid (2.41 bard)

Sanitization & Sterilization

Autoclave	250 °F (121 °C), 30 min, 5+ cycles
Chemical Sanitization	Performed using industry standard concentrations of hydrogen peroxide, peracetic acid, sodium hypochlorite and other selected chemicals.

Endotoxins

The levels of bacterial endotoxins in aqueous extracts from BCWPS Micro Capsule filters are below current USP limits as specified for water for injection.

Extractables

BCWPS Micro Capsule filters typically exhibit low levels of non-volatile residues.

TOC and Conductivity

BCWPS Micro Capsule filters conform with TOC standards of USP <643> and the water conductivity standards of USP <645> after an appropriate flush with purified water.

Toxicity Compliance

Materials used to construct BCWPS Micro Capsule filters are non-toxic and meet the requirements for the MEM Elution Cytotoxicity Test and the requirements for Biological Reactivity Tests in the current version of the United States Pharmacopeia (USP) for Class VI - 121 °C Plastics.

Non-Fiber Releasing

BCWPS Micro Capsule filters comply with Title 21 CFR sections 210.3 (b)(6) and 211.72, for non-fiber releasing filters.

FDA Compliance

Materials meet the requirements listed by the FDA as appropriate for use in articles intended for repeated food contact as specified in Title 21 CFR sections 174.5, 177.1500, 177.1520, 177.1630, 177.2440, and 177.2600 as applicable.

Single Layer Filtration Area (Nominal)

Area	0.569 ft ²
	529 cm ²


Double Layer Filtration Area (Nominal)

Area	0.45 ft ²
	418 cm ²

Construction Materials

Filtration Media	Single or Dual Layer High Capacity PES Membrane with Polyester Support
Media Support	Polypropylene
End Caps, Center Core, Outer Support Cage, Micro Capsule Housing	Polypropylene
Sealing Method	Thermal Bonding

Flow Rates for BCWPS Micro Capsules by Pore Size



Flow rates for Micro Capsule filters are per filter. The test fluid is water at ambient temperature. Flows are tested using a Micro capsule filter with ½" Sanitary inlet and outlet ports. Rates will vary based on end configuration of the Micro capsule.

BCWPS Micro Capsule Filters Ordering Information

All Critical Process filters are configurable to meet customer specifications.
Fill in the corresponding codes in the boxes below to build your Part Number.

To consult with one of our technical team members, request a quote or place an order:
call (603) 880-4420 or [contact us here](#).

Micro Capsule Filters

MIC BCWPS

Prefilter Pore Size Code Blank = no prefilter -20 = 0.2 µm -50 = 0.5 µm -80 = 0.8 µm 1-0 = 1.0 µm		Final Layer Pore Size Code -10 = 0.1 µm -20 = 0.2 µm -50 = 0.5 µm -80 = 0.8 µm 1-0 = 1.0 µm		Inlet 1 = 1/8" Hose Barb 2 = 1/4" Hose Barb 3 = 1/2" Hose Barb 4 = Luer Lock 5 = 1/2" Sanitary* 6 = 1/4" MNPT		Inlet Vent Port 1 = 1/8" Hose Barb 2 = 1/4" Hose Barb 3 = 1/2" Hose Barb 4 = Luer Lock 6 = 1/4" MNPT 7 = Side Bleed Valve		Outlet 1 = 1/8" Hose Barb 2 = 1/4" Hose Barb 3 = 1/2" Hose Barb 3B = 1/2" Hose Barb with Filling Bell 4 = Luer Lock 5 = 1/2" Sanitary* 6 = 1/4" MNPT		Outlet Vent Port 1 = 1/8" Hose Barb 2 = 1/4" Hose Barb 3 = 1/2" Hose Barb 4 = Luer Lock 6 = 1/4" MNPT 7 = Side Bleed Valve		Side Vent O-Ring** B = Buna E = EP S = Silicone V = Viton (or FKM) K = FFKM
		Pre-Sterilized or Not S = Pre-Sterilized N = Not Sterilized										

*When choosing the Sanitary Inlet/Outlet, the Luer Lock or bleed valve option is required for the Vent Port
** O-Ring is only available on Bleed Valve



One Chestnut Street
Nashua, NH 03060
603.880.4420
FAX: 603.880.4536

CriticalProcess.com

The information contained herein is subject to change without notice. The Critical Process Filtration logo is a trademark of Critical Process Filtration, Inc. Viton is a trademark of DuPont Performance Elastomers L.L.C.
© 2025 Critical Process Filtration, Inc. • All Rights Reserved

Data Sheet BCWPS Micro DS Rev -