



3X $\phi 10\text{mm}$ ($\phi 39.4''$) EQUALLY SPACED
 DN A $\phi 350\text{mm}$ ($\phi 13.78''$) BC

DIMENSIONS & SANITARY PIPE CONNECTIONS						
"A"	NO.	NAME	PORT SIZE	QTY	MATERIAL	STYLE
10" - 24.0"	N1	INLET	2"	1	316L	CLAMP
20" - 34.0"	N2	OUTLET	2"	1	316L	CLAMP
30" - 44.0"	N3	GAUGE PORT	1"	1	316L	CLAMP
40" - 54.0"	N4/N5	DRAIN PORTS	0.5"	1	316L	CLAMP
	N6	VENT PORT	0.5"	1	316L	CLAMP

HOUSING SPECIFICATIONS	
DESIGNED PRESSURE (PSIG)	195
OPERATING PRESSURE (PSIG)	150
OPERATING TEMP (°C)	110

- TECHNICAL PARAMETERS**
1. THIS APPARATUS MEETS GMP STANDARD
 2. ALL WELD-SEAMS ARE ARGON WELDED
 3. MAXIMUM CORROSION RESISTANCE
 4. SURFACE FINISH: 400# MACHINE POLISH-PASSIVATE-ELECTROPOLISH ALL SURFACES.
 5. MATERIAL: 316L

UNLESS OTHERWISE SPECIFIED		MATERIALS		DATE	
All dimensions are in inches. Remove all burrs and break sharp edges. Tolerances: Decimals		SURFACE TREATMENT		DWN BY	
.XX ± .01 .XXX ± .005 Angles ± 1°		The drawing and the data contained hereon are the property of Critical Process, Nashua, NH. U.S.A. It is to be used in manufacturing and for reproduction of this drawing is expressly prohibited without Critical Process' consent.		BMD/DN 3/08	
DO NOT SCALE PRINT		ACAD FILE # CSH3.DWG		DATE	
		LIST OF MATERIAL		CHK BY	
		NUMBER		APRD BY	
				DATE	
				APRD BY	
				DATE	
				RELEASE	
				DATE	
				TITLE	
				SALES DRAWING	
				3 ROUND	
				SANITARY T-STYLE	
				SITE	
				D	
				DWG NO	
				CSH3	
				REV	
				-	

CRITICAL PROCESS FILTRATION, INC
 NASHUA, NH
 WWW.CRITICALPROCESS.COM
 PH: 603-880-4420