

Use this guide to find the end fittings that will fit your connections. The codes to use for your inlet and outlet fittings when ordering your capsule filter are shown in the list below. Reference our individual product data sheets to find detailed ordering information.



1/4" Female NPT



1/4" Male NPT



3/8" Female NPT



1/2" Female NPT



1/2" Male NPT



Sanitary\*



Hose Barb\*\*

## End Fitting Codes

### Inlet

A = 1/4" Female NPT  
B = 1/4" Male NPT  
C = 3/8" Female NPT  
D = 1/2" Female NPT  
E = 1/2" Male NPT  
F = Sanitary\*  
G = Hose Barb\*\*

### Outlet

A = 1/4" Female NPT  
B = 1/4" Male NPT  
C = 3/8" Female NPT  
D = 1/2" Female NPT  
E = 1/2" Male NPT  
F = Sanitary\*  
G = Hose Barb\*\*

\* Fits 1" or 1 1/2" sanitary gaskets/clamps

\*\*Hose Barb Diameter Ranges

	Minimum	Maximum
<b>Outer Diameters</b>	11/32" (8.6mm)	9/16" (14.0mm)
<b>Inner Diameters</b>	5/32" (4.0mm)	13/32" (10.5mm)

Visit our [website](#) or [contact us](#) for more information and to access datasheets on all of our products.