

Capsule Filter End Fitting Guide



Use this guide to find the end fittings that will fit your connections. The codes to use for your inlet and outlet fittings when ordering your capsule filter are shown in the list below. Reference our individual product data sheets to find detailed ordering information.



1/4" Female NPT



1/4" Male NPT



3/8" Female NPT



1/2" Female NPT



1/2" Male NPT



1"-1 1/2" Sanitary*



Hose Barb**

*Fits 1" or 1 1/2" sanitary gaskets/clamps

**Tapered Hose Barb Diameter Ranges:

	Minimum	Maximum
Outer Diameters	11/32" (8.6mm)	9/16" (14.0mm)
Inner Diameters	5/32" (4.0mm)	13/32" (10.55mm)

End Fitting Codes

Inlet

A = 1/4" Female NPT
 B = 1/4" Male NPT
 C = 3/8" Female NPT
 D = 1/2" Female NPT
 E = 1/2" Male NPT
 F = 1"-1 1/2" Sanitary*
 G = Hose Barb**

Outlet

A = 1/4" Female NPT
 B = 1/4" Male NPT
 C = 3/8" Female NPT
 D = 1/2" Female NPT
 E = 1/2" Male NPT
 F = 1"-1 1/2" Sanitary*
 G = Hose Barb**