



Mini-Capsule Filter End Fitting Guide

Use this guide to find the end fittings that will fit your connections. On the Mini-Capsule's datasheet, fill in the appropriate end fitting code listed below when creating your part number. Photos represent outlet configurations with Luer Lock vent ports.



1/8" Hose Barb



1/4" Hose Barb



1/2" Hose Barb



1/2" Hose Barb with Filling Bell



Luer Lock



1/2" Sanitary*



1/4" MNPT

*When choosing the sanitary inlet/outlet, the Luer Lock option is required for the Vent Port

End Fitting Codes

Inlet

- 1 = 1/8" Hose Barb
- 2 = 1/4" Hose Barb
- 3 = 1/2" Hose Barb
- 4 = Luer Lock
- 5 = 1/2" Sanitary*
- 6 = 1/4" MNPT

Inlet Vent Port

- 1 = 1/8" Hose Barb
- 2 = 1/4" Hose Barb
- 3 = 1/2" Hose Barb
- 4 = Luer Lock
- 6 = 1/4" MNPT

Outlet

- 1 = 1/8" Hose Barb
- 2 = 1/4" Hose Barb
- 3 = 1/2" Hose Barb
- 3B = 1/2" Hose Barb with filling bell
- 4 = Luer Lock
- 5 = 1/2" Sanitary*
- 6 = 1/4" MNPT

Outlet Vent Port

- 1 = 1/8" Hose Barb
- 2 = 1/4" Hose Barb
- 3 = 1/2" Hose Barb
- 4 = Luer Lock
- 6 = 1/4" MNPT



One Chestnut Street
 Nashua, NH 03060
 603.880.4420
 FAX: 603.880.4536
 CriticalProcess.com

The information contained herein is subject to change without notice. The Critical Process Filtration logo is a trademark of Critical Process Filtration, Inc.
 © 2021 Critical Process Filtration, Inc. • All Rights Reserved

Mini-Capsule End Fitting Guide Rev-