



DIMENSIONS & SANITARY PIPE CONNECTIONS						
"A"	NO.	NAME	PORT SIZE	QTY	MATERIAL	STYLE
10" - 29.6"	N1	INLET	2"	1	316L	CLAMP
20" - 39.6"	N2	OUTLET	2"	1	316L	CLAMP
30" - 49.6"	N3	GAUGE PORT	1"	1	316L	CLAMP
40" - 59.6"	N4/N5	DRAIN PORTS	0.5"	1	316L	CLAMP
	N6	VENT PORT	0.5"	1	316L	CLAMP

HOUSING SPECIFICATIONS	
DESIGNED PRESSURE (PSIG)	195
OPERATING PRESSURE (PSIG)	150
OPERATING TEMP (°C)	110

- TECHNICAL PARAMETERS**
1. THIS APPARATUS MEETS GMP STANDARD
 2. ALL WELD-SEAMS ARE ARGON WELDED
 3. MAXIMUM CORROSION RESISTANCE
 4. SURFACE FINISH: 400# MACHINE POLISH - PASSIVATE-ELECTROPOLISH ALL SURFACES.
 5. MATERIAL: 316L

UNLESS OTHERWISE SPECIFIED		MATERIALS		DATE	
All dimensions are in inches. Remove all burrs and break sharp edges. Tolerances: Decimals		SURFACE TREATMENT		3/08	
XX ± .01 XXX ± .005 Angles ± 1°		The drawing and the data contained hereon are the property of Critical Process, Nashua, NH. U.S.A. It is to be used in manufacturing and for reproduction of the drawing is expressly prohibited without Critical Process consent.		DWG NO	
Surface Roughness		APPRO BY		DATE	
DO NOT SCALE PRINT		RELEASE		DATE	
		TITLE		SALES DRAWING	
		NASHUA, NH		6 ROUND	
		WWW.CRITICALPROCESS.COM		SANITARY T-STYLE	
		PH: 603-880-4420		CSH6	
		CRITICAL PROCESS FILTRATION, INC		REV	
		ACAD FILE # CSH6.DWG		NUMBER	
		LIST OF MATERIAL			
		ITEM		QTY	
		MATERIAL			