


DIMENSIONS & SANITARY PIPE CONNECTIONS						
"A"	NO.	NAME	PORT SIZE	QTY	MATERIAL	STYLE
5" - 12.6"	N1	INLET	1.5"	1	316L	CLAMP
10" - 17.6"	N2	OUTLET	1.5"	1	316L	CLAMP
20" - 27.6"	N3	GAUGE PORT	0.5"	1	316L	CLAMP
30" - 37.6"	N4	DRAIN PORT	0.5"	1	316L	CLAMP
40" - 47.6"						

HOUSING SPECIFICATIONS	
DESIGNED PRESSURE (PSIG)	195
OPERATING PRESSURE (PSIG)	150
OPERATING TEMP (°C)	110



- TECHNICAL PARAMETERS**
1. THIS APPARATUS MEETS GMP STANDARD
 2. ALL WELD-SEAMS ARE ARGON WELDED
 3. MAXIMUM CORROSION RESISTANCE
 4. SURFACE FINISH: 400# MACHINE POLISH-PASSIVATE-ELECTROPOLISH ALL SURFACES.
 5. MATERIAL: 316L

ITEM	QTY	LIST OF MATERIAL	NUMBER
ACAD FILE # CSH1.DWG			
UNLESS OTHERWISE SPECIFIED		MATERIALS	
All dimensions are in inches. All burns and holes are 1/8" deep. Tolerances:		SURFACE TREATMENT	
Decimals: .XX ± .01		APRD BY DATE	
Angles: ± 1°		APRD BY DATE	
Surface Roughness: 43		APRD BY DATE	
DO NOT SCALE PRINT		RELEASE DATE	
<p>The design and the data contained hereon are confidential and the exclusive property of Critical Process, Nashua, NH. No part of this drawing may be reproduced or transmitted in any form or by any means, electronic, mechanical, photocopying, recording, or by any information storage and retrieval system, without the prior written permission of Critical Process.</p>			
CRITICAL PROCESS FILTRATION, INC		TITLE	
NASHUA, NH		SALES DRAWING	
WWW.CRITICALPROCESS.COM		SANITARY INLINE	
PH: 603-880-4420		VENT HOUSING	
D		CSH1	
REV		-	