

Capsule Filter End Fitting Guide

Use this guide to find the end fittings that will fit your connections. The codes to use for your inlet and outlet fittings when ordering your capsule filter are shown in the list below. Reference our individual product data sheets to find detailed ordering information.



1/4" Female NPT



1/4" Male NPT



3/8" Female NPT



1/2" Female NPT



1/2" Male NPT



1" Sanitary*



Stepped Hose Barb**



1 1/2" Sanitary with Side Vent



1/2" Single Stepped Barb with Side Vent



3/4" Single Step Barb with Side Vent



Compression (JACO®)



6mm Quick Disconnect

*Fits 1" or 1/2" Sanitary Gaskets/Clamps

**Tapered Hose Barb Diameter Ranges

	Minimum	Maximum
Outer Diameter	11/32" (8.6mm)	9/16" (14.0mm)
Inner Diameter	5/32" (4.0mm)	13/32" (10.55mm)

End Fitting Codes and Options

Inlet

- A = 1/4" FNPT
- B = 1/4" MNPT
- C = 3/8" FNPT
- D = 1/2" FNPT
- E = 1/2" MNPT
- F = 1" Sanitary
- G = Stepped Hose Barb
- H = 1 1/2" Sanitary with Side Vent
- I = 1/2" Single Stepped Barb with Side Vent
- J = 3/4" Single Step Barb with Side Vent
- Y = Compression (JACO®)
- Z = 6mm Quick Disconnect

Outlet

- A = 1/4" FNPT
- B = 1/4" MNPT
- C = 3/8" FNPT
- D = 1/2" FNPT
- E = 1/2" MNPT
- F = 1" Sanitary
- G = Stepped Hose Barb
- H = 1 1/2" Sanitary with Side Vent
- I = 1/2" Single Stepped Barb with Side Vent
- IB = 1/2" Single Stepped Barb with Side Vent & Filling Bell
- J = 3/4" Single Step Barb with Side Vent
- Y = Compression (JACO®)
- Z = 6mm Quick Disconnect