

Micro Capsule Filter End Fitting Guide

Use this guide to find the end fittings that will fit your connections. On the Micro Capsule's datasheet, fill in the appropriate end fitting code listed below when creating your part number. Photos represent outlet configurations with Luer Lock vent ports.



1/8" Hose Barb



1/4" Hose Barb



1/2" Hose Barb



Luer Lock



1/2" Sanitary*



1/4" MNPT

*When choosing the sanitary inlet/outlet, the Luer Lock option is required for the Vent Port

End Fitting Codes and Options

Inlet

- 1 = 1/8" Hose Barb
- 2 = 1/4" Hose Barb
- 3 = 1/2" Hose Barb
- 4 = Luer Lock
- 5 = 1/2" Sanitary*
- 6 = 1/4" MNPT

Inlet Vent Port

- 1 = 1/8" Hose Barb
- 2 = 1/4" Hose Barb
- 3 = 1/2" Hose Barb
- 4 = Luer Lock
- 6 = 1/4" MNPT

Outlet

- 1 = 1/8" Hose Barb
- 2 = 1/4" Hose Barb
- 3 = 1/2" Hose Barb
- 3B = 1/2" Hose Barb with Filling Bell
- 4 = Luer Lock
- 5 = 1/2" Sanitary*
- 6 = 1/4" MNPT

Outlet Vent Port

- 1 = 1/8" Hose Barb
- 2 = 1/4" Hose Barb
- 3 = 1/2" Hose Barb
- 4 = Luer Lock
- 6 = 1/4" MNPT